



Essederm<sup>®</sup>

Medium

**Essederm<sup>®</sup> Medium**

Is a Sterile restorable cross-linked pyrogen free hyaluronic acid gel of non-animal origin which is used for: Filling medium and deep cutaneous depressions such as wrinkles, folds and other defects, restoring or creating volume in facial contours (cheekbones, chin) and lip contours

- \* 1.0 ml cross link hyaluronic acid
- \* 29G disposable glass syringe
- \* It is sterile pyrogen free, colorless, isotonic gels, transparent, odorless, viscostic, no particle and 100% none animal base
- \* Cross linked NaHA content 11.5 to 12.5 mg per ml
- \* Not cross linked NaHA content up to 1.2 mg per ml
- \* Cross linking grade: up to 90 %
- \* DVS Content: < 0.02 mg per ml
- \* Storage and condition: 2 to 25 C (35 to 77 F)

Essederm<sup>®</sup> Medium

Family of products

Rev.01

Essederm®

1. Worry Lines, Essederm **Medium**
2. Glabellar Lines, Essederm **Medium**
4. Augmentation of Chickbones, Essederm **Medium**
5. Cheek Augmentation / Cheek Wrinkles, Essederm **Medium**
6. Nasolabial Folds, Essederm **Medium**
8. Vermilion Border, Essederm **Medium**
9. Oral Commissurer / Smile Lines, Essederm **Medium**
10. Lip Contour, Essederm **Medium**
12. Chin Augmentation, Essederm **Medium**

