



Essederm[®]

Subline

Essederm[®] Subline

Is a Sterile restorable cross-linked pyrogen free hyaluronic acid gel of non-animal origin which is used for: deep skin deformations and depressions, increasing volume and contours of the face (chin, cheek bones), rhinoplasty and scar correction

- * 1.0 ml cross link hyaluronic acid
- * 26G disposable glass syringe
- * It is sterile pyrogen free, colorless, isotonic gels, transparent, odorless, viscous, no particle and 100% none animal base
- * Cross linked NaHA content 26 to 28 mg per ml
- * Not cross linked NaHA content up to 1.2 mg per ml
- * Cross linking grade: up to 90 %
- * DVS Content: < 0.02 mg per ml
- * Storage and condition: 2 to 25 C (35 to 77 F)

Essederm[®] Subline

Family of products

Rev.01

Essederm®



4. Augmentation of Chickbones, Essederm Subline
5. Cheek Augmentation / Cheek Wrinkles, Essederm Subline
6. Nasolabial Folds, Essederm Subline

12. Chin Augmentation, Essederm Subline

