



Essederm[®]

Thin

Essederm[®] Thin

Is a Sterile restorable cross-linked pyrogen free hyaluronic acid gel of non-animal origin which is used for:
Filling superficial wrinkles (perioral, periorbital...),
Augmentation of lip and injection in lip contour

- * 1.0 ml cross link hyaluronic acid
- * 30 G disposable glass syringe
- * It is sterile pyrogen free, colorless, isotonic gels, transparent, odorless, viscous, no particle and 100% none animal base
- * Cross linked NaHA content 5.6 to 6.5 mg per ml
- * Not cross linked NaHA content up to 1.2 mg per ml
- * Cross linking grade: up to 90 %
- * DVS Content: < 0.02 mg per ml
- * Storage and condition: 2 to 25 C (35 to 77 F)

Essederm[®] Thin

Family of products
Rev.01

Essederm®



- 1. Worry Lines, Essederm **Thin**
- 2. Glabellar Lines, Essederm **Thin**
- 3. Periorbital Lines, Essederm **Thin**

7. Perioral Lines, Essederm **Thin**

11. Lip Augmentation, Essederm **Thin**

