

Essederm®




1. Worry Lines, Essederm **Thin, Medium**
2. Glabellar Lines, Essederm **Thin, Medium**
3. Periorbital Lines, Essederm **Thin**
4. Augmentation of Chickbones, Essederm **Medium, Deep, Subline**
5. Cheek Augmentation / Cheek Wrinkles, Essederm **Medium, Deep, Subline**
6. Nasolabial Folds, Essederm **Medium, Deep, Subline**
7. Perioral Lines, Essederm **Thin**
8. Vermilion Border, Essederm **Medium**
9. Oral Commissurer / Smile Lines, Essederm **Medium, Deep**
10. Lip Contour, Essederm **Medium**
11. Lip Augmentation, Essederm **Thin**
12. Chin Augmentation, Essederm **Medium, Deep, Subline**



Essederm®

Family of products



The highest concentration of cross link up to 90% (Up to 4 times higher than others)
Each package contains 1 x 1.0 ml syringes

	Concentration	Needle size	Indications	Skin level Injection Area
<i>Essederm</i> Thin	6 mg/ml	30 G	Superficial Wrinkles Crow's feet Red part of the lips Perioral lines	Superficial Dermis A
<i>Essederm</i> Medium	12 mg/ml	29 G	Medium and Deep Wrinkles In all areas Lip contours	Middle and Deep Dermis B & C
<i>Essederm</i> Deep	18 mg/ml	27 G	Deep Wrinkles Contour of the face Volume creation	Deep Dermis or Subcutis C & D
<i>Essederm</i> Subline	27 mg/ml	26 G	Specially Contour Of the face Very Deep Wrinkles Volume creation	Subcutis and Epiperiostal E
Injection technique	based on the indications: Retrotracing, Criss Cross, Multiple Point and Fan Technique			

Rev. 01